

Healthpac Server Shutdown Procedure

Prior to shutting down the Healthpac server, it is important to shutdown the Healthpac database first. The following is a step by step procedure for this process. This can be done from any workstation with access to H2000 or from the server.

```
Enter your user id: DM
```

```
terminal name: MM-MON
product name: IBM VIDEO
terminal width: 79      printer width: 80
                depth: 24          depth: 59
lineskip:      0
lf delay:      1
ff delay:      1
back space:    8
```

```
[1301] The Speller is off.
master dictionary:DM
```

```

                                     3
<<< 11:59:18  --                    D                    --    16 Feb
2005  >>>
<<<
>>>
<<< Copyright (c) 1973-1999 Pick Systems.  All rights
reserved. >>>
<<< This work is property of, and embodies trade secrets
and    >>>
<<< confidential information proprietary to Pick Systems.
It    >>>
<<< may not be reproduced, copied, used, disclosed,
transferred >>>
<<< adapted or modified without the express written
approval of >>>
<<< Pick Systems.
>>>
```

```
:SHUTDOWN
```

Do you wish to continue (Y/N/O)?**Y**

No host defined.

Terminating all print jobs...

Wrapping up process on line 0

Wrapping up process on line 2

Wrapping up process on line 3

Wrapping up process on line 4

Wrapping up process on line 5

Wrapping up process on line 6

Wrapping up process on line 7

Wrapping up process on line 8

Wrapping up process on line 9

Wrapping up process on line 10

Wrapping up process on line 11

Wrapping up process on line 12

Wrapping up process on line 13

Wrapping up process on line 14

Wrapping up process on line 15

Wrapping up process on line 16

Wrapping up process on line 17

Wrapping up process on line 18

Wrapping up process on line 19

Wrapping up process on line 20

Wrapping up process on line 21

Wrapping up process on line 22

Wrapping up process on line 23

Wrapping up process on line 24

Wrapping up process on line 25

Wrapping up process on line 26

Wrapping up process on line 27

Wrapping up process on line 28

Wrapping up process on line 29

Wrapping up process on line 30

Wrapping up process on line 31

Wrapping up process on line 32

Wrapping up process on line 37

Wrapping up current process on line 1

Pause for wrapup processing to complete...

Flushing memory, boot when disk is quiescent...

Once you are disconnected from Accuterm, you may shutdown or reboot the Healthpac server.

Please call 800-831-9419 option1 for questions.